# INSTRUMENTS

# PCE-DHM 5 DIGITAL HANDHELD MICROSCOPE Deutsch/English/日本語/中文/Español/Français

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Thanks for choosing our microscope. Please read this manual carefully and keep it in a safe place. If you hand the product to others, please be sure to attach this manual. Please heed all warnings and follow all instructions in the manual.

### Microscope Kit Introduction

PCE-DHM 5 is a new concept of digital microscope that can be used as a handheld or tabletop magnification tool, to provide bright and clear images, it offers a 6-LED illumination system with adjustable bright white output. To offer flexibility in viewing visual data, the microscope features a tiltable 2" color TFT LCD screen, which can be tilted 75" the most convenient viewing position. and Long-term observation is comfortable for your cervical spine. That can be a great portable companion in your home, outdoor, lab, school, work, or in the field. が疲れません。

#### Microscope Kit Checklist (This microscope kit includes the following parts) Microscope and Accessories Use Description

(1) The digital Microscope – The microscope holds samples of very small things on clear glass slides. It shines a bright light from its bed plate LED lamp through the small sample,and then lenses make the sample look very big.  $(\widehat{2})$  Type-C data cable - used to transmit data, charge the microscope, and charge the observation platform.

(3) Cleaning cloth----only used to wipe some large areas of dirt or water droplets on the microscope body.

(4) Carrying case and Straps – This is designed for safe storage of the microscope, shock resistance for transportation, and portability for the

(5) Display - This is where you view and set up and perform microscope shooting functions. The image projected by the acquisition objective lens is includes the parameter setting of the display and the brightness adjustment of the microscope illumination. 8 - button 9 Power button for detailed operation, please refer to the microscope display buttons function and interface description. (10) Focus Ring - Turn this ring very slowly to focus on the specimen you are viewing so you can see it clearly. (11) Objective - It makes the observed object play a role in visual magnification. (12) Transparent plastic cover - mainly for you to clearly see the current state of the specimen when you observe the specimen (the height of the cover is just the best focusing distance of the microscope) (13) LEDs Ring Light - Illuminates the observed object to make it clearly visible.

(14) Output Ports - Includes Type-C USB port and TF card slot, allowing you to charge the microscope and store image/video data. (please cover the 15 Microscope charging indicator 16 TF card slot (17) Type-C charging port (18) Reset hole

 Long press the (9) on the monitor to light up the screen. 1.2. Please long press the (7) to see if the light is turned on, and short press the (7) (8) to see if the zoom changes. Press and hold the (8) to see if 1. Aim the microscope at the object to be observed, make sure the object is centered on the circular light-transmitting cover ()), and at the edge of 2-2. Adjust the inclination of the display screen to suit your viewing angle  $\widehat{(2)}$ \*The display can be only adjusted at the tilt angle (about 75°) for viewing, please do not rotate it horizontally, otherwise, it may damage the hinge of 2-3. Please gently rotate the focus ring of the microscope, and carefully observe the screen changes. If it is clear, you can take pictures or videos  $\mathfrak{D}/\mathfrak{Q}$ . For more operations, please refer to "Setup Instructions"

Microscope display buttons function and interface description Press and hold the power button (9) for 2 seconds, the startup screen will appear. In any mode, press and hold the power button (9) for 2 seconds to turn off the device. on the device, the default is photo mode. Short press the power buttoh @ to switch the modes @ <math>2 / 2 / 2.

Short press the + button to zoom in, the 7 to zoom out, and the 9 to take pictures/videos. Press and hold the 7 to turn on the light/brightness+, (8) brightness-/turn off the light. Long press the (6) to enter the settings menu for this (29/30)/(3).

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