Deutsch

English

# Bedienungsanleitung User Manual

PCE Smart Measure Environmental

User manuals in various languages (français, italiano, español, português, nederlands, türk, polski, русский, 中文) can be found by using our product search on: [www.pce-instruments.com](http://www.pce-instruments.com)

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## English Contents

<b>1</b>	<b>Safety notes</b> .....	<b>12</b>
<b>2</b>	<b>System requirements</b> .....	<b>12</b>
2.1	iOS.....	12
<b>3</b>	<b>Installation</b> .....	<b>12</b>
<b>4</b>	<b>Description of the user interface</b> .....	<b>13</b>
<b>5</b>	<b>Meaning of the individual icons of the main menu</b> .....	<b>14</b>
<b>6</b>	<b>Operation</b> .....	<b>16</b>
6.1	First use of the app .....	16
6.2	Connect to the meter.....	16
6.3	Disconnect from the meter .....	17
<b>7</b>	<b>Make a measurement</b> .....	<b>17</b>
7.1	PCE-MSL 2BT and PCE-HWA 20BT .....	17
7.2	PCE-PMI 1BT .....	17
7.3	Representation of the readings.....	18
7.4	Language.....	18
<b>8</b>	<b>Share measurement data</b> .....	<b>19</b>
8.1	Save measurement data .....	20
8.2	Export series of measurements .....	22

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## 1 Safety notes

Acquaint yourself with the operation of the measuring device before using it in combination with the app. For this purpose, use the instruction manual which comes with your PCE Instruments meter. All safety notes from the meter's manual have to be observed as well when the device is used with this app.

The device and the app may only be used by qualified personnel. Damage or injuries caused by non-observance of the manuals are excluded from our liability and not covered by our warranty.

We do not assume liability for printing errors or any other mistakes in this manual.

We expressly point to our general warranty terms which can be found in our general terms of business.

If you have any questions, please contact PCE Instruments. The contact details can be found at the end of this manual.

## 2 System requirements

### 2.1 Android

- Android version 9 Pie (API 28) or higher
- Bluetooth interface (version 4.2)
- Screen size 5.71 inches
- A minimum resolution of 1520x720 pixels
- Processor: ARM Cortex-A53, 2000 Mhz, 4 cores
- 2 GB RAM recommended

### 2.2 iOS


- Current iOS version
- Bluetooth interface (version 4.2)
- Screen size 5.8 inches
- 2GB RAM

## 3 Installation

Please download the app from the Google Play Store or Apple App Store and install it on your smartphone. Then check and grant the access permissions for the location and memory.

## 4 Description of the user interface

The main window is composed of two sections. In the upper toolbar, there is a menu button on the left which opens a side menu. The side menu contains menu items for software information and company information. More detailed information about these menu items will follow below.




Below the toolbar, there are six buttons, each representing a functionality.

## 5 Meaning of the individual icons of the main menu

Main menu	
	Open window for Bluetooth connections
	Open window for measurements
	Open window for sharing
	Open window for data
	Open window for help
	Open window for export
	Open side menu
"Connection" window	
	Device disconnected
	Device connected
Search devices	Search for Bluetooth devices
Disconnect	Disconnect device
	Back to main menu
"Measure" window	
	Initiate new measurement
	Back to main menu
Save	Save measured value (PCE-PMI 1BT)
"Share" window	
	Back to main menu
	Share selection via messenger

 Measurement 8/16/2021 09:49:50 <input checked="" type="checkbox"/>	Select a measurement
<b>"Data" window</b>	
	Back to main menu
 Measurement 8/16/2021 09:49:50  	Open window for measurement data
	Rename data record
	Delete data record
<b>"Export" window</b>	
 Measurement 8/16/2021 09:49:50 <input checked="" type="checkbox"/>	Select a measurement
	Export selection
	Back to main menu
<b>Side menu</b>	
  About	Open window for software and company information


## 6 Operation

### 6.1 First use of the app

Before the meter can work with the app, make sure that Bluetooth is activated. Furthermore, the access permissions for location and memory must be granted. These are used exclusively to search for Bluetooth devices in the direct vicinity and to save the measurement data to the smartphone as PDF files. The permissions can be granted, depending on the device, via **Settings -> Apps -> PCE Smart Measure Environmental -> Permissions**. After the access permissions have been granted, the app can be used to its full extent.

### 6.2 Connect to the meter

After starting the app, a Bluetooth connection to the meter can be established. It must be ensured that the smartphone and the meter are not more than 5... 10 m apart. To establish the connection, navigate to the corresponding menu by tapping on "**Connection**" in the main menu. Bluetooth connections can be managed under "**Connection**". To establish a connection, a suitable meter must be searched for. This can be initiated by tapping on the "**Search devices**" button. The devices found are then displayed in a list, see figure.

Once the meter has been found (e. g. PCE-HWA 20BT), a connection can be established by tapping on the meter in the list. After the connection has been successfully established, the user immediately receives a confirmation on the display that the devices are connected and is forwarded directly to the main menu.

### 6.3 Disconnect from the meter

Tap on the  icon to disconnect the meter. Closing the app while the Bluetooth connection is active disconnects the devices. After use and before closing the mobile application, the devices must always be disconnected as described above. This ensures correct functionality of the device.

## 7 Make a measurement

The menu item "**Measure**", which can be reached from the main menu, offers the possibility to perform a measurement. Before making the first measurement, ensure that there is an active Bluetooth connection between the smartphone and the meter. Since the meters have different features, a measurement is made differently depending on the meter you are using.

### 7.1 PCE-MSL 2BT and PCE-HWA 20BT


With the meters PCE-MSL 2BT and PCE-HWA 20BT, the measured values are saved when viewing mode is opened. To initiate a new measurement, click on the  icon in the menu bar. You can then name the measurement and confirm with OK.

### 7.2 PCE-PMI 1BT

When using the moisture meter PCE-PMI 1BT, you will find the  button in the app. Tap on this button to save a value to the database. As it is done with the two meters mentioned above, click on the  icon in the menu bar to initiate a new measurement.

### 7.3 Representation of the readings

Depending on the meter, different measured values are shown with their units.


English

### 7.4 Language

Depending on the language set in your smartphone, the language and format will be displayed either in English or German.

## 8 Share measurement data


All performed measurements can be shared via all common messengers as a PDF file. To do so, call up the menu item "**Share**" in the main menu. The window contains a list of all measurements that have been made. The desired measurements can be marked accordingly by placing check marks.

After tapping on the  button, you can select the desired export format. An options dialogue opens to select the desired messenger. After the selection, the chosen measurements are attached as PDF or CSV documents and can be sent.

## 8.1 Save measurement data

All performed measurements are saved to a database in the memory of the smartphone. These can be viewed, edited or deleted at any time. To do this, call up the menu item "Data" in the main menu. The view contains all measurement data named according to the time stamp of the measurement. As the mobile application can be used with several meters, you must select the meter the measurements of which you want to view.

Each measurement can be renamed or deleted via the associated options. Furthermore, it is possible to view the measured values by tapping on the desired measurement.

The screenshot shows a mobile application interface titled "Data Details". It features a table with four columns: "Time", "Temp. [°C] ↑", "RH [%]", and "Speed [m/s]". The table contains 13 rows of data, showing a sequence of measurements over time. The temperature starts at 0.0 °C and then jumps to 22.0 °C. Relative humidity (RH) is 0.0% for the first two rows and then increases to 60.8% and 60.9% for the remaining rows. The speed is consistently 0.0 m/s or 0.1 m/s.


Time	Temp. [°C] ↑	RH [%]	Speed [m/s]
12:06:42	0.0 °C	0.0 %	0.0 m/s
12:06:42	22.0 °C	0.0 %	0.0 m/s
12:06:42	22.0 °C	60.8 %	0.1 m/s
12:06:43	22.0 °C	60.8 %	0.1 m/s
12:06:43	22.0 °C	60.9 %	0.1 m/s
12:06:43	22.0 °C	60.9 %	0.1 m/s
12:06:43	22.0 °C	60.9 %	0.1 m/s
12:06:44	22.0 °C	60.9 %	0.1 m/s
12:06:44	22.0 °C	60.9 %	0.1 m/s
12:06:44	22.0 °C	60.9 %	0.1 m/s
12:06:44	22.0 °C	60.9 %	0.1 m/s
12:06:45	22.0 °C	60.8 %	0.1 m/s

In the overview of the saved measurement data, the individual values are shown in a chart. Here, you can sort the readings by their value by tapping on the relevant column header.

## 8.2 Export series of measurements

Series of measurements can be exported via the menu item "**Export**" in the main menu.

In the "**Export**" window, individual series of measurements can be selected and then exported by tapping on the  button. A progress bar is displayed while the data are being processed for export. After processing, the series of measurements are saved to the "**PCE**" folder of the smartphone file system as a PDF or CSV file.

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Specifications are subject to change without notice.

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PCE Americas Inc.
PCE Teknik Cihazları Ltd.Şti.