

Sample Cutter MPC-3



PCE Americas Inc.
711 Commerce Way
Suite 8
Jupiter
FL-33458
USA
From outside US: +1
Tel: (561) 320-9162
Fax: (561) 320-9176
info@pce-americas.com

PCE Instruments UK Ltd.
Unit 11
Southpoint Business Park
Ensign way
Hampshire / Southampton
United Kingdom, SO31 4RF
From outside UK +44
Tel: (0) 2380 98703 0
Fax: (0) 2380 98703 9
info@pce-instruments.co.uk

www.pce-instruments.com/english
www.pce-instruments.com

Delivery content

- 1x rubber pad

Spare parts:

- MPC-EK3 spare blades
(1 package = 40 pieces)
- MPC-UL-N rubber pads for
sample cutter (10 pieces)

Sample cutter MPC-3

Mechanical sample cutter for the preparation of paper, fleece and textile samples to determine the surface weight.

- sample cutter for cutting thicknesses of up to 3 mm
- sample surface of 100 square centimetres (suitable for the paper scales).
- accuracy of the cut-out 5%

The sample cutters are used to cut materials like paper, cardboard, fleece, fabric, knitted fabric, foil, felt, all types of textiles, etc.

Operation of the sample cutter:

Put the sample material on a cutting board and place the sample cutter on the sample material. Hold the device with one hand and release the safety mechanism below the rotary handle (pull, then turn) and turn the upper rotary handle by 90°.

