

## Technical Crane Scale PCE-PS 300MLS

### **Crane Scale PCE-PS 300MLS (M class III)** **Crane Scale up to 300 kg / easy-to-read touch-screen / tare function / powered by mains component or accumulator / shackle included / RS-232 port and optional software**

The PCE-PS 300MLS Crane Scale has a wide weighing range up to 300 kg and can be integrated in your present bariatric system. The PCE-PS 300MLS Crane Scale has an easy-to-use display by means of two keys to on/off the scale, zero set or tare. The PCE-PS 300MLS Crane Scale is powered by a mains component or an internal accumulator which are included in the delivery. The internal accumulator offers up to 80 hours of continuous operation. The PCE-PS 300MLS Crane Scale is the ideal solution to weigh people who cannot move themselves but it is necessary to control the patient's weight and it is not possible to use a conventional scale. Many bariatric systems have connections and fixing devices to place a bariatric scale. The PCE-PS 300MLS verifiable Crane Scale can be integrated in many bariatric systems. Crane Scales are very demanding in nursing residences and hospitals. The PCE-PS 300MLS Crane Scale is send with the verification done. The first verification is valid for two years and must be recalibrated annually. If you have any questions regarding our bariatric scales, please contact us and our engineers and technicians will help you to find the Crane Scale that best suits your needs.

- Shackle included in the delivery
- Easy-to-read display
- Easy-to-use keyboard
- Resolution of 100 g
- Weighing range: 300 kg
- HOLD function
- Tare along the entire range
- RS-232 port / optional software
- 230 V mains adapter
- Accumulator for 80 h of continuous operation
- Verifiable according to the commercial class M III

#### Technical Specifications

Weighing range	300 kg
Readability	100 g
Verification value	100 g
Minimum weight	2 kg
Verification class	III
Tare range	Along the entire range
Weighing units	kg
Response time	2 s

Display	LCD (with backlight)
Mainsadapter	230 V / 50 Hz
Accumulator	12 V / 100 mA
Dimensions	210 x 130 x 190 mm
Operatingtemperature	+10 ... +40 °C
Weight	approx. 2,6 kg

**Delivery contents**

1 x PCE-PS 300MLS Crane Scale, 1 x accumulator, 1 x mains component, user's manual