

## Hoja de Datos

Mezclar / Agitadores Magnéticos



### color squid IKAMAG® white

The improved small magnetic stirrers now in new designs.

new: Digital speed display (LED)

new: Electronically controlled motor for more capacity

new: Higher speed range from 0 - 2.500 rpm

new: Max. stirring quantity 1 l

- Outstanding chemical resistance due to glass top and synthetic bottom made of TPC-ET

- Recyclable materials

Accesorios: IKAFLON® 15 Varillas magnéticas de agitación, IKAFLON® 30 Varillas magnéticas de agitación, IKAFLON® 10 Varillas magnéticas de agitación, IKAFLON® 25 Varillas magnéticas de agitación, TRIKA® 25 Varillas magnéticas, RSE Extractor de varillas agitadoras, IKAFLON® 20 Varillas magnéticas de agitación

#### Información Técnica

Número de puestos de agitación	1
Max. cantidad a agitar por cada posición (H2O) [l]	1
Max.cantidad a agitar (H2O) [l]	1
Potencia del consumo del motor [W]	3
Potencia del suministro del motor [W]	2
Indicador de velocidad	LED
Rango de velocidad [rpm]	0 - 2500
Longitud max. varilla magnética [mm]	30
Dimensiones de la superficie de la placa [mm]	Ø 115
Dimensiones (An x Al x Pr) [mm]	145 x 45 x 160
Peso [kg]	0.54
Temperatura ambiental permitida [°C]	5 - 40
Humedad relativa permitida [%]	80
Clase de protección de acuerdo al DIN EN 60529	IP 54
Interfaz de RS 232	no
Potencia análoga	no
Voltaje [V]	100 - 240
Frecuencia [Hz]	50/60
Consumo de energía [W]	4
Consumo de energía [mA]	300
<b>No. ref.</b>	<b>3671000</b>