

Hoja de Datos

Mezclar / Agitadores Magnéticos



big squid IKAMAG® white

The improved magnetic stirrers now in new designs.

new: Digital speed display (LED)

new: Electronically controlled motor for more capacity

new: Higher speed range from 0 - 2.500 rpm

- Outstanding chemical resistance due to glass top and synthetic bottom made of TPC-ET

- Recyclable materials

Accesorios: IKAFLON® 15 Varillas magnéticas de agitación, IKAFLON® 30 Varillas magnéticas de agitación, IKAFLON® 10 Varillas magnéticas de agitación, IKAFLON® 25 Varillas magnéticas de agitación, TRIKA® 25 Varillas magnéticas, RSE Extractor de varillas agitadoras, IKAFLON® 20 Varillas magnéticas de agitación

Información Técnica

Número de puestos de agitación	1
Max. cantidad a agitar por cada posición (H2O) [l]	1.5
Max.cantidad a agitar (H2O) [l]	1
Potencia del consumo del motor [W]	3
Potencia del suministro del motor [W]	2
Indicador de velocidad	LED
Rango de velocidad [rpm]	0 - 2500
Longitud max. varilla magnética [mm]	30
Control de velocidad	50 pasos de RPM
Material de la superficie de la placa	vidrio
Dimensiones de la superficie de la placa [mm]	Ø 160
Dimensiones (An x Al x Pr) [mm]	180 x 40 x 195
Peso [kg]	0.7
Temperatura ambiental permitida [°C]	5 - 40
Humedad relativa permitida [%]	80
Clase de protección de acuerdo al DIN EN 60529	IP 54
Interfaz de RS 232	no
Potencia análoga	no
Frecuencia [Hz]	50/60
Consumo de energía [W]	4
Consumo de energía [mA]	300
No. ref.	3672000