

PCE-BSH series RS-232 OUTPUT FORMAT

Modus :KEY / STB / ON / OFF

Baud Rate :1200、2400、4800、9600

Data Bit :8/7/7

Parity : None/ Even/Odd

Stop Bit :1

Code :ASCII

Data Format :

W=WEIGHT M=MESSAGE

1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
HEAD		DATA											UNIT				CR		
W/M	:	+/-											g					CR	LF

Example:

100.00g : W:+ 100.00g

500.00ct : W:+ 500.00ct

Err-0 : M:over

data cable

