



# Paper Scale Round Sample Cutter PCE-FSC 10



## **Paper scale round sample cutter for 10 cm<sup>2</sup> samples, robust and compact, four blades, up to 5 mm**

The paper scale round sample cutter impresses with its compact dimensions and the robust steel housing. With the four internal blades of the round sample cutter, you can cut out a 10 cm<sup>2</sup> sample (fleece, paper, leather, textile, cork, ...) within seconds, a quarter turn is sufficient, and you can measure its weight per unit area.

- ▶ Robust and compact
- ▶ Four internal stainless steel blades
- ▶ Up to 0-5 mm / 0-0.19" material thickness
- ▶ 10 cm<sup>2</sup> / 1.55 in<sup>2</sup> sample size
- ▶ With mechanical lock
- ▶ Suitable for many materials

Subject to change

## Specifications

Sample thickness	max. 0-5 mm / 0-0.19"
Sample cut-out	<b>10 cm<sup>2</sup> / 1.55 in<sup>2</sup></b>
Housing	cast steel
Cutting blades	stainless steel
Dimensions	126 mm / 4.96", Ø 84 mm / 3.30"
Weight	ca. 500 g / 17.6 oz

## More information

More product info



Similar products



Subject to change