

Photometer PCE-TUM 50







Photometer PCE-TUM 50

Turbidity meter for drinking and waste water as well as food / Colored graphic display /Measurement of NTU, EBC, ASBC / According to ISO 7027 (90°) / AC Mains operation (220V AC power adapter) / USB interface

The photometer PCE-TUM 50 is a portable device for measuring turbidity according to ISO 7027. This allows the user to quickly and easily determine important water parameters with the turbidity meter and to check the quality of drinking water and waste water. The calibration of the turbidity meter PCE-TUM 50 is very easy due to the calibration standards supplied. The calibration can be carried out with 7 preset points. In addition, user-specific calibration points can also be set. The turbidity meter is capable of measuring turbidity of liquids up to 2000 NTU.

Furthermore, the measuring unit can be switched to EBC (European Brewery Convention) or ASBC (American Society of Brewing Chemists). The measured values are displayed on the large color display. The turbidity meter is mainly used in the laboratory and is therefore operated with a plug-in power supply. Via the rear USB interface, the turbidity meter can be connected to a computer and the memory can be cleared.

What is the EBC or European Brewery Convention? (See above in the picture gallery)As EBC, the color or the color intensity of beer and wort is described in Europe. This set value describes how much light is absorbed by beer with a certain original wort content. The turbidity meter can determine this value and thus helps in judging the type of beer.

- ► High precision
- ▶ 3 measuring units
- ▶ Turbidity calibration up to 7 points
- ► Color display
- ► Easy handling
- ▶ USB interface



Specifications

Measurement method ISO 7027 (90°)

Measuring ranges 0 ... 2000 NTU

0 ... 500 EBC 0 ... 9999 ASBC

Resolution < 100 NTU: 0.01

< 999 NTU: 0.1 < 2000 NTU: 1

Accuracy < 500 NTU: ± 2% of measured value

< 2000 NTU: ± 3% of measured value

Calibration < 0.02 NTU

10 NTU
200 NTU
500 NTU
1000 NTU
1500 NTU
2000 NTU

Or user-specific

Light source 850 nm IR LED Scattered light influence < 0.02 NTU

Dimensions of cuvettes 60 x 25 mm / 2.4 x 1 in

Automatic shutdown 2 hours after the device was not used

Storage 200 memory slots

Interface USB

Power supply 220 VAC/50 Hz, 12 VDC

Current consumption About 300-mA

Operating temperature 0 ... 60°C / 32 ... 140°F

Dimensions 250 x 177 x 96 mm / 9.8 x 7 x 3.8 in

Weight About 1.2 kg / 2.6 lbs

More information

Manual

More product info



Similar products

